


**1. Part Number:** A856-100

**2. Photograph:**



**3. Product Information:**

Description

Stainless Container with 6" Top 100L

Material of Construction:

Body: 316L stainless steel  
 Lid: 316L stainless steel  
 Clamp: 316L stainless steel (may be 304 stainless steel)  
 Gasket: Silicone (conforms to FDA 21 CFR 177.2600)

Method of Construction:

Fully welded, the weld is ground and polished

Surface Finish:

Better than 0.5 microns Ra

Dimensions:

A = 167mm  
 B = 150mm  
 C = 670mm  
 D = 465mm  
 E = 500mm  
 Wall Thickness = 2mm

Nominal Volume:

100L

Graduation Marks:

None

Pressure Rating:

None

Nominal Weight (body, lid, clamp & gasket):

20.9 Kg

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin

Storage Conditions:

Dry, ambient conditions

Note: The neck of the container is fitted with a 6" diameter ferrule which conforms to BS4825-3

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Feb-22