


**1. Part Number:** A856-75

**2. Photograph:**



**3. Product Information:**

|   |   |
|---|---|
| Description                                 | Stainless Container with 6" Top 75L   |
| Material of Construction:                   | Body: 316L stainless steel<br>Lid: 316L stainless steel<br>Clamp: 316L stainless steel (may be 304 stainless steel)<br>Gasket: Silicone (conforms to FDA 21 CFR 177.2600) |
| Method of Construction:                     | Fully welded, the weld is ground and polished   |
| Surface Finish:                             | Better than 0.5 microns Ra  |
| Dimensions:                                 | A = 167mm<br>B = 150mm<br>C = 650mm<br>D = 465mm<br>E = 450mm<br>Wall Thickness = 2mm   |
| Nominal Volume:                             | 75L   |
| Graduation Marks:                           | None  |
| Pressure Rating:                            | None  |
| Nominal Weight (body, lid, clamp & gasket): | 17.7 Kg   |
| BSE/TSE statement:                          | All polishing compounds used in the manufacture of this product are of vegetable origin   |
| Storage Conditions:                         | Dry, ambient conditions   |

Note: The neck of the container is fitted with a 6" diameter ferrule which conforms to BS4825-3

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Feb-22



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.  
 Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com