




1. Part Number: 1030CT all sizes

2. Product Information:



Part No:	Nominal Volume	Weight
1030CT-10	10ml	299g
1030CT-20	20ml	221g
1030CT-25	25ml	182g

Description

Tips for use with 25mm diameter PowderThief Bodies

Materials of Construction:

316L stainless steel

Finish:

<1 microns Ra

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin

Recommended Storage Conditions:

Dry and ambient temperature

3. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25