

1. Part Number: 1095N-150

2. Photo:



### 3. Product Information:

Description	V-Spatula 150mm x 20mm
Method of Construction:	Crevice free, all welds ground and polished
Materials:	316L stainless steel
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	233mm
Blade Length:	150mm
Maximum Blade Width:	20mm
Nominal Weight:	150g
Recommended Storage Conditions:	Dry and ambient temperature
First Use:	This sampler should be thoroughly cleaned before being used for the first time

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Nov-19