

## **Specification Sheet**

Specification Number: 1210K-2000
Revision: 03
Date Issued: 4-Oct-19

**1. Part Number:** 1210K-2000

2. Photo:




## 3. Product Information:

Description Jumbo Slot Sampler 2000mm

Material of Construction: 316L stainless steel

Outer Diameter of Outer Tube: 40 mm
Inner Diameter of Inner Tube: 32 mm
Overall Length: 2060 mm
Slot Length: 120 mm
Slot Width: 20 mm
Number of Slots

Slot Positions:



External Polish: <1 Micron Ra

Internal Polish: None

BSE/TSE statement: All polishing compounds used in the manufacture of this product are

of vegetable origin

Nominal Weight of Sampler: 6.5 Kg

Recommended Storage Conditions: Dry and ambient temperature

4. Approval: Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Oct-19



4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com