

1. Part Number: 1213A-01

2. Photo:



3. Product Information:

Description	Powder Max 1000mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Internal Polish:	None
Material of Construction:	316L stainless steel
Overall Sampler Length:	1040mm
Outer Diameter of Tube:	38mm
Inner Diameter of Tube:	34mm
Maximum Diameter of Sleeve:	52mm
Slot Width:	30mm
Slot Length:	30mm
Sampling Slot Position:	



Nominal Weight:	2 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-24