

**1. Part Number:** 1640A-1500

**2. Photograph:**



### 3. Product Information:

Description:	Point Sampler 1500mm
Material of Construction:	316L stainless steel
Method of Construction:	Welds ground and polished
Finish	- External Bright Polish - Internal Machine finish
Nominal Sample Volume:	100 ml
Overall Length:	1730 mm
Max. Diameter:	42 mm (NB. The scraper will add an extra 10mm)
Nominal Weight of Sampler:	2.1 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this item are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-20