

1. Part Number: 1900A-10

2. Photo:



3. Product Information:

Description	Measuring Scoop 10ml
Nominal Volume:	10ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Outer Scoop Diameter:	25mm
Inner Scoop Diameter:	22mm
Scoop Depth:	30mm
Overall Length:	162mm
Nominal Weight:	90g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 15-Nov-24