

**1. Part Number:** 5005A-500

**2. Photo:**



### 3. Product Information:

|                             |                                                                                         |
|-----------------------------|-----------------------------------------------------------------------------------------|
| Part No.                    | 5005A-500                                                                               |
| Description                 | Stainless Steel Dipper – 500ml                                                          |
| Material of Construction:   | 316L stainless steel                                                                    |
| Method of Construction:     | Fully Welded Construction                                                               |
| Surface Finish:             | Cup better than 0.5 microns Ra<br>Handle better than 1.0 microns Ra                     |
| Nominal Cup Volume:         | 500ml                                                                                   |
| Brimful Cup Volume:         | 500ml                                                                                   |
| Overall Length:             | 600mm                                                                                   |
| Max. Width (cup and handle) | 86mm                                                                                    |
| Cup Diameter:               | 76mm                                                                                    |
| Cup Height:                 | 120mm                                                                                   |
| Nominal Weight:             | 730g                                                                                    |
| BSE/TSE statement:          | All polishing compounds used in the manufacture of this product are of vegetable origin |

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Nov-19