

**1. Part Number:** 5020L-04S**2. Product Information:**

Description: Polyethylene Tubing 4m – Sterile  
Tubing Length: 4m (approx.)  
Quantity of Tubes per Box 15  
Material of Construction: Low Density Polyethylene (LDPE)  
EC Compliance 1935/2004, 2023/2006, 10/2011, 202/2014  
US Compliance FDA 21 CFR 177.1520  
Individually bagged Yes  
Sterile Yes (Gamma Irradiation 25-45kGy)  
Overall Diameter: 8mm  
Internal Diameter: 6mm  
Pressure Rating: None  
ATEX: This product has not been ATEX tested  
Box Weight: 1.8 Kg  
Recommended Storage Conditions: Dry and ambient temperature

**3. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-24