

**1. Part Numbers:** 5551S-1500

**2. Photograph:**



**3. Product Information:**

Description	Fluid Point Sampler 1500mm
Materials of Construction:	<ul style="list-style-type: none"> <li>- 316L stainless steel</li> <li>- Silicone O-ring (FDA acceptable CFR 177.2600)</li> </ul>
Diameter of Outer Tube:	25mm
Diameter of hole in tip:	2mm
Max. Length:	3050mm
Min Length:	1595mm
Nominal Weight:	2.8kg
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Dec-24