

**1. Part Number:** 8005A-02

**2. Photograph:**



**3. Product Information:**

Description:	Sealing Sleeve 40mm
Use:	For use with 40mm diameter Pocket Samplers
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"><li>• Conforms to FDA CFR 177.1520</li><li>• Conforms to EU Regulations 10/2011</li><li>• Conforms to EC Regulations 1935/2004</li><li>• BSE/TSE free</li></ul>
Method of Construction:	Single piece
Max. Sampler Diameter:	40mm
Sleeve Length:	160mm
Nominal Weight:	24g
Number of Sleeves per Pack:	3
Recommended Storage Conditions:	Dry and ambient temperature

Note: The Pocket Sampler shown in the photograph is for illustrative purposes only. The sampler must be purchased separately to the Sealing Sleeves

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 7-Oct-20