

Specification Sheet

Specification Number: 8060H-510S
Revision: 12
Date Issued: 21-Nov-22

1. Part Number: 8060H-510S

2. Photograph:



3. Product Information:

Description: SteriWare PowderThief - 500mm – 10ml Tip - Sterile

Use: Single use, disposable powder sampler

Material of Construction: HDPE (High Density Polyethylene) - virgin

• Conforms to FDA CFR 177.1520

Conforms to EU Regulations 10/2011
Conforms to EC Regulations 1935/2004

BSE/TSE free

Method of Construction: Ultrasonic welding

Nominal Length: 511mm

Diameter of Outer Tube: 21mm

Nominal Tip Size: 10ml

Width of Handle: 103mm

Length 'L': 63mm

Nominal Weight of Sampler: 93g

Moulding Environment: Class 100,000 Cleanroom
Assembly & Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (minimum dose 25 kGy)

Number of Samplers per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture
ATEX: This device has not been ATEX tested





Specification Sheet

Specification Number: 8060H-510S Revision: 12 Date Issued: 21-Nov-22

4. Change Control: Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the

SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin: Made and packed in the UK.

6. Repeated Use: This is a single-use, disposable sampler. Repeated use will lead to

failure.

7. Approval: Signed:

Name: Kuldeep Flora
Position: Quality Manager

Date: 21-Nov-22