

**1. Part Number:** 8060H-520

**2. Photograph:**



### 3. Product Information:

|   |  |
|---|--|
| Description:                            | SteriWare PowderThief - 500mm – 20ml Tip   |
| Use:                                    | Single use, disposable powder sampler  |
| Material of Construction:               | HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul> |
| Method of Construction:                 | Ultrasonic welding   |
| Nominal Length:                         | 511mm  |
| Diameter of Outer Tube:                 | 21mm   |
| Nominal Tip Size:                       | 20ml   |
| Width of Handle:                        | 103mm  |
| Length 'L':                             | 121mm  |
| Nominal Weight of Sampler:              | 84g  |
| Moulding Environment:                   | Class 100,000 Cleanroom  |
| Assembly & Primary Packing Environment: | Class 10,000 Cleanroom   |
| Individually Bagged?                    | Yes (heat sealed PE/PET bag)   |
| Method of Sterilisation:                | N/A  |
| Number of Samplers per Box:             | 20   |
| BSE/TSE Free:                           | Yes  |
| Recommended Storage Conditions:         | Dry and ambient temperature  |
| Shelf Life:                             | 5 years from the date of manufacture   |
| ATEX:                                   | This device has not been ATEX tested   |



**SAMPLING  
SYSTEMS**

# Specification Sheet

**Specification Number:** 8060H-520  
**Revision:** 10  
**Date Issued:** 22-Nov-22



**ISO  
9001 : 2015  
REGISTERED**  
Cert No. 14128455

**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK.  
Tel: +44(0)1675 466992 Email: [info@sampling.com](mailto:info@sampling.com) Web: [www.sampling.com](http://www.sampling.com)

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



**Name:** Kuldeep Flora  
**Position:** Quality Manager  
**Date:** 22-Nov-22