

Specification Sheet

Specification Number: 8081P-3000S Revision: 07

Date Issued: 23-Nov-22

1. Part Number: 8081P-3000S

2. Photograph:



3. Product Information:

Description: SteriWare Tray 3000ml - sterile

Use: Single use, disposable tray

Material of Construction: Polypropylene (PP) - rigid

Conforms to FDA CFR 177.1520

Conforms to EU Regulation 10/2011

• Conforms to EC Regulation 1935/2004

BSE/TSE free

Method of Construction: Single Piece, injection moulded

Overall Dimensions: 318 x 267 x 55mm

Nominal Volume: 3000ml Nominal Weight: 180g

Moulding Environment:

Class 100,000 cleanroom.

Primary Packing Environment:

Class 10,000 cleanroom.

Yes (heat sealed PE bag)

Method of Sterilisation:

Gamma Irradiation (25-45kGy)

Number of Trays per Box: 30

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 2 years from the date of manufacture
ATEX: This device has not been ATEX tested



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com



Specification Sheet

Specification Number: 8081P-3000S Revision: 07 Date Issued: 23-Nov-22

4. Change Control: Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin: Made and packed in the UK.

6. Repeated Use: This is a single-use, disposable product. Repeated use will lead to

failure.

7. Approval: Signed:

Name: Kuldeep Flora

Quality Manager

Date: 23-Nov-22

Position: