

1. Part Number: 8098A-25S

2. Photograph:



3. Product Information:

Description:	SteriWare® Weighing Boat 25ml - Sterile
Use:	Single use, disposable weighing boat
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> • FDA acceptable CFR 177.1640 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	25ml
Overall Length:	93mm
Maximum Width:	59mm
Depth:	13mm
Nominal Weight:	2g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45 kGy)
Quantity per Box:	30

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Nov-24