

1. Part Number: 8098A-50S

2. Photograph:



3. Product Information:

| | |
|------------------------------|--|
| Description: | SteriWare® Weighing Boat 50ml - Sterile |
| Use: | Single use, disposable weighing boat |
| Material of Construction: | High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> • FDA acceptable CFR 177.1640 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Nominal Volume: | 50ml |
| Overall Length: | 115mm |
| Maximum Width: | 72mm |
| Depth: | 18mm |
| Nominal Weight: | 4g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Primary Packing Environment: | Class 10,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25 to 45 kGy) |

Quantity per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 12-Nov-24