

**1. Part Number:** 8098A-50

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Weighing Boat 50ml
Use:	Single use, disposable weighing boat
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1640</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC regulation 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Volume:	50ml
Overall Length:	115mm
Maximum Width:	72mm
Depth:	18mm
Nominal Weight:	4g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A

Quantity per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 12-Nov-24