

**1. Part Numbers:** 8126H-500A

**2. Photograph:**



### 3. Product Information:

|                                 |                                                                                                           |
|---------------------------------|-----------------------------------------------------------------------------------------------------------|
| Description                     | Wide Mouth HDPE Bottle 500ml Amber                                                                        |
| Materials of Construction:      | Bottle HDPE (conforms to FDA 177.1520 & USP Class VI)<br>Cap PP (conforms to FDA 177.1520 & USP Class VI) |
| Manufacturing Environment:      | Controlled environment within factory                                                                     |
| Nominal Volume:                 | 500ml                                                                                                     |
| Brimful Volume:                 | 550ml                                                                                                     |
| Nominal Height:                 | 169mm                                                                                                     |
| Diameter:                       | 73mm                                                                                                      |
| Internal Neck Diameter:         | 43mm                                                                                                      |
| Weight:                         | 68g                                                                                                       |
| Pack Size:                      | 48 bottles                                                                                                |
| Method of Sterilisation:        | N/A (these bottles are non-sterile)                                                                       |
| BSE/TSE Free:                   | Yes                                                                                                       |
| Supplied with lid screwed on:   | Yes                                                                                                       |
| Recommended Storage Conditions: | Dry and ambient temperature                                                                               |

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24