

1. Part Numbers: 8136H-500

2. Photograph:



3. Product Information:

Description	Wide Mouth HDPE Bottle 500ml
Materials of Construction:	Bottle = HDPE Cap = PP
	All polymers conform to: • FDA 21 CFR 177.1520 • USP Class VI <87> & <88> • EU 10/2011
Manufacturing Environment:	ISO Class 8 Cleanroom (Class 100,000)
Nominal Volume:	500ml
Brimful Volume:	537ml
Nominal Height:	165mm
Diameter:	75mm
Internal Neck Diameter:	43mm
Weight:	72g
Pack Size:	48 bottles per pack
Method of Sterilisation:	N/A (these bottles are non-sterile)
BSE/TSE Free:	Yes
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Endotoxin Free:	Yes - tested to ISO 11737-3 & USP <85>
Particulate Free:	Yes - tested to USP <788>



Phthalate Free: Yes

Leakproof: Yes – No leaks observed when a sample bottle was filled with water, inverted with a pressure of 715mmHg for 15 minutes

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Aug-26