

**1. Part Number:** 8185P-500

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® GMP Tub 500ml
Use:	Wide mouth bottle with screw on cap - for use with powders
Material of Construction:	Polypropylene (PP) – virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulation 10/2011</li> <li>• Conforms to EC Regulation 1935/2004</li> <li>• BSE/TSE free</li> </ul>
	White Colouring (1%)
Method of Construction:	Single piece, injection moulded
Nominal Volume:	500ml
Brimful Volume:	540ml
Height:	106mm
Lid Diameter:	94mm
Nominal Weight:	83g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	No
Method of Sterilisation:	N/A
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

ATEX:

This product has not been ATEX tested

#### **4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

#### **5. Country of Origin:**

Made and packed in the UK.

#### **6. Approval:**

Signed:



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 23-Nov-22