

1. Part Number: 8209H-05S

2. Photograph:



3. Product Information:

| | |
|---------------------------------|--|
| Description: | SteriWare® Long Handled Spoon 5ml (Straight) - Sterile |
| Use: | Sampling and dispensing spoon |
| Material of Construction: | High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Bowl Volume: | 5ml |
| Handle Length: | 285mm |
| Bowl Length: | 46mm |
| Bowl Width: | 30mm |
| Nominal Weight: | 17g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Primary Packing Environment: | Class 10,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (min dose 25kGy) |
| Quantity per Box: | 100 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | This product has not been ATEX tested |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Mar-26