


1. Part Number: 8284A-150

2. Photograph:



3. Product Information:

Description:	SteriWare® Cupped Spatula – 150mm
Use:	Single use, disposable spatula
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> • FDA acceptable CFR 177.1640 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Blade Shape:	U
Blade Style:	Closed Tip
Nominal Volume (flat filled):	21ml
Nominal Length:	240mm
Blade Length:	155mm
Blade Width:	18mm
Nominal Weight:	20g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture

ATEX:

4. Change Control:

This item has not been ATEX tested

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com**5. Country of Origin:**

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Aug-24