

# Product Specification

Specification Number

8404H-1000S

Revision

02

Part No.	8404H-1000S		
Description	Wide Mouth HDPE Bottle 1000ml - Sterile		
Materials of Construction:	Bottle	HDPE	(conforms to FDA 177.1520 & USP Class VI)
	Cap	PP	(conforms to FDA 177.1520 & USP Class VI)



Manufacturing Environment:	Controlled environment within factory
Packing Environment:	Class 100,000 Cleanroom
Nominal Volume:	1000ml
Brimful Volume:	1075ml
Nominal Height:	201mm
Diameter:	92mm
Internal Neck Diameter:	52mm
Weight:	100g
Individually Bagged:	Yes (in a PE/PET laminate film bag)
Pack Size:	30 bottle per box
Method of Sterilisation:	Gamma Irradiation (25-45 kGy)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes

Specification Approved By:

Signature:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-May-19