

Specification Sheet

Specification Number: 8432P-1000S Revision: 03

Date Issued: 25-Nov-22

1. Part Number: 8432P-1000S

2. Photograph:



3. Product Information:

Description: SteriWare® Graduated Jug 1000ml – Sterile

Use: Single use, disposable jug

Material of Construction: Polypropylene (PP) – virgin

• Conforms to FDA CFR 177.1520

• Conforms to EU Regulations 10/2011

• Conforms to EC Regulations 1935/2004

BSE/TSE free

Method of Construction: Single piece, injection moulded

Nominal Volume: 1000ml
Brimful Volume: 1270ml
Top Diameter (excl. spout): 120mm
Max Width (inc. spout and handle) 170mm
Base Diameter: 93mm
Height: 170mm
Nominal Weight: 123g

Graduation Marks: To 1000ml (every 10ml)

Moulding Environment: Class 100,000 Cleanroom

Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25 – 45 kGy)



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com



Specification Sheet

Specification Number: 8432P-1000S Revision: 03 Date Issued: 25-Nov-22

Quantity per Box: 18

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture
ATEX: This device has not been ATEX tested

4. Change Control: Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the

SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin: Made and packed in the UK.

6. Repeated Use: This is a single-use, disposable sampler. Repeated use will lead to

failure.

7. Approval: Signed:

Name: Kuldeep Flora
Position: Quality Manager

Date: 25-Nov-22