

**1. Part Number:** 8432P-500

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Graduated Jug 500ml
Use:	Single use, disposable jug
Material of Construction:	Polypropylene (PP) – virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Volume:	500ml
Brimful Volume:	635ml
Top Diameter (excl. spout):	90mm
Max Width (inc. spout and handle)	135mm
Base Diameter:	75mm
Height:	130mm
Nominal Weight:	58g
Graduation Marks:	To 500ml (every 10ml)
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)

Method of Sterilisation: N/A  
Quantity per Box: 35  
Recommended Storage Conditions: Dry and ambient temperature  
Shelf Life: 5 years from the date of manufacture  
ATEX: This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.  
To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 25-Nov-22