


**1. Part Number:** 8601S-02

**2. Photograph:**



### 3. Product Information:

**Description:**

Stainless Sample Splitter – 6mm slots

**Material of Construction:**

Sample Splitter: Stainless steel  
Collecting Bins: Stainless steel  
Frame: Painted mild steel

**Method of Construction:**

Spot welded

**Finish:**

Unpolished, welds left as laid

**Number of slots:**

16

**Slot Width:**

6mm

**Overall Height:**

335mm

**Bin Dimensions:**

163 x 140 x 84mm

**Number of Bins Supplied:**

2

**Nominal Weight:**

1.90 kg (Splitter & 2 Bins)

**Recommended Storage Conditions:**

Dry and ambient temperature

### 4. Approval:

**Signed:**



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 27-Oct-20