

1. Part Number: 8601S-10

2. Photograph:



3. Product Information:

Description:	Stainless Sample Splitter – 10mm slots	
Material of Construction:	Sample Splitter:	Stainless steel
	Collecting Bins:	Stainless steel
	Frame:	Painted mild steel
Method of Construction:	Spot welded	
Finish:	Unpolished, welds left as laid	
Number of slots:	16	
Slot Width:	10mm	
Overall Height:	348mm	
Bin Dimensions:	273 x 137 x 84mm	
Number of Bins Supplied:	2	
Nominal Weight:	2.2kg (Splitter & 2 Bins)	
Recommended Storage Conditions:	Dry and ambient temperature	

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 27-Oct-20