


1. Part Number: 8601S-20

2. Photograph:



3. Product Information:

| | | |
|---------------------------------|--|--------------------|
| Description: | Stainless Sample Splitter – 20mm slots | |
| Material of Construction: | Sample Splitter: | Stainless steel |
| | Collecting Bins: | Stainless steel |
| | Frame: | Painted mild steel |
| Method of Construction: | Spot welded | |
| Finish: | Unpolished, welds left as laid | |
| Number of slots: | 12 | |
| Slot Width: | 20mm | |
| Overall Height: | 390mm | |
| Bin Dimensions: | 343 x 137 x 84mm | |
| Number of Bins Supplied: | 2 | |
| Nominal Weight: | 3.20kg (Splitter & 2 Bins) | |
| Recommended Storage Conditions: | Dry and ambient temperature | |

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20