


1. Part Numbers: 8644E-100

2. Photograph:



3. Product Information:

Description:	Cleaning Coupon EPDM – 100mm
Material of Construction:	EPDM – Polymer conforms to: • FDA 21 CFR 177.2600 • USP Class VI • BSE/TSE free
Finish:	n/a
Dimensions:	100 x 100mm
Nominal Weight:	40g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name:	Kuldeep Flora
Position:	Quality Manager
Date:	23-Aug-24