


 SAMPLING SYSTEMS	Product Specification	Specification Number	8702A-1
		Revision	01

Part No. 8702A-1

Description Thermometer

Material of Construction: Stainless steel probe
ABS plastic controls with LCD display



Diameter of Probe: 3mm

Probe Length: 130mm

Overall Length: 246mm

Nominal Weight of Thermometer: 45g

Range -39.9 to 149.9°C


Accuracy -39.9 to 124.9°C (+/- 0.5°C)

Display LCD

Battery Life 3000 hours
(Auto power off after 10 minutes)

Replacement Battery: 3V Lithium Coin Cell Type CR2032

Calibration Certificate: Yes

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 24-Jun-14
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	