

 SAMPLING SYSTEMS	Sampler Specification	Specification Number	8703A-25
		Revision	01

Part No. 8703A-25

Description Titrating Automatic Burette 25ml

Product Contact Materials: Polyethylene
Borosilicate Glass
Latex



Method of Construction: Parts push together, they can be dis-assembled for cleaning


Overall Height: 61mm

Capacity of Burette: 25ml

Graduation Marks on Burette: Every 0.10ml
(Accuracy +/- 0.045ml, Category B ISO 385)

Volume of Reservoir Bottle: 1000ml

Nominal Weight: 194g

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 24-Jun-14
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	