

 SAMPLING SYSTEMS	<h1 style="text-align: center;">Sampler Specification</h1>	Specification Number	8703A-50
		Revision	01

Part No. 8703A-50

Description Titrating Automatic Burette 50ml

Product Contact Materials: Polyethylene
Borosilicate Glass
Latex



Method of Construction: Parts push together, they can be dis-assembled for cleaning


Overall Height: 870mm

Capacity of Burette: 50ml

Graduation Marks on Burette: Every 0.10ml
(Accuracy +/- 0.075ml, Category B ISO 385)

Volume of Reservoir Bottle: 1000ml

Nominal Weight: 225g

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 24-Jun-14
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	