

## **Specification Sheet**

Specification Number: 8785B-75S Revision: 03

25-Nov-22

Date Issued:

**1. Part Number:** 8785B-75S

### 2. Photograph:



#### 3. Product Information:

Description: SteriWare® Scraper 75mm Blue - Sterile

Use: Single use, disposable scraper

Material of Construction: Polypropylene – blue - virgin

• Conforms to FDA CFR 177.1520

• Conforms to EU Regulations 10/2011

• Conforms to EC Regulations 1935/2004

• BSE/TSE free

Method of Construction: Single piece, injection moulded

Overall Length: 190mm
Blade Length: 75mm
Max. Blade Width: 75mm
Nominal Weight: 36g

Moulding Environment: Class 100,000 Cleanroom
Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (Minimum 25 kGy)

Quantity per Box: 100

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture
ATEX: This device has not been ATEX tested



### Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com



# **Specification Sheet**

Specification Number: 8785B-75S
Revision: 03
Date Issued: 25-Nov-22

4. Change Control: Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the

SteriWare range please email your contact details to

change.control@sampling.com

**5. Country of Origin:** Made and packed in the UK.

**6. Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to

failure.

7. Approval: Signed:

Name: Kuldeep Flora

**Quality Manager** 

Date: 25-Nov-22

Position: