

**1. Part Number:** A123-100

**2. Photograph:**



### 3. Product Information:

Description	GMP Pallet Knife 100mm
Material:	Stainless Steel - grade 420
Surface Finish:	- Better than 0.5 microns Ra - Crevice Free
Overall Length:	188mm
Blade Length:	100mm
Max. Blade Width:	20mm
Nominal Weight:	64g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature
Cleaning:	This product should be thoroughly cleaned before first use. Use the same detergents and cleaning protocols that are used for other small items of stainless steel processing equipment.

### 4. Approval:

Signed:



**Name:** Kuldeep Flora  
**Position:** Quality Manager  
**Date:** 17-Dec-24