

1. Part Number: A123-200

2. Photograph:



3. Product Information:

Description	GMP Pallet Knife 200mm
Material:	Stainless Steel - grade 420
Surface Finish:	- Better than 0.5 microns Ra - Crevice Free
Overall Length:	325mm
Blade Length:	210mm
Max. Blade Width:	31mm
Nominal Weight:	150g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature
Cleaning:	This product should be thoroughly cleaned before first use. Use the same detergents and cleaning protocols that are used for other small items of stainless steel processing equipment.

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 17-Dec-24