

1. Part Number: A123-300

2. Photograph:



3. Product Information:

| | |
|---------------------------------|--|
| Description | GMP Pallet Knife 300mm |
| Material: | Stainless Steel - grade 420 |
| Surface Finish: | - Better than 0.5 microns Ra - Crevice Free |
| Overall Length: | 442mm |
| Blade Length: | 300mm |
| Max. Blade Width: | 51mm |
| Nominal Weight: | 290g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Cleaning: | This product should be thoroughly cleaned before first use. Use the same detergents and cleaning protocols that are used for other small items of stainless steel processing equipment. |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 17-Dec-24