


**1. Part Number:** A123-350

**2. Photograph:**



**3. Product Information:**

Description	GMP Pallet Knife 350mm
Material:	Stainless Steel - grade 420
Surface Finish:	<ul style="list-style-type: none"> <li>- Better than 0.5 microns Ra</li> <li>- Crevice Free</li> </ul>
Overall Length:	490mm
Blade Length:	350mm
Max. Blade Width:	51mm
Nominal Weight:	290g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature
Cleaning:	<p>This product should be thoroughly cleaned before first use.</p> <p>Use the same detergents and cleaning protocols that are used for other small items of stainless steel processing equipment.</p>

**4. Approval:**

Signed:



**Name:** Kuldeep Flora  
**Position:** Quality Manager  
**Date:** 17-Dec-24