

**1. Part Numbers:** A125-1000

**2. Photograph:**



**3. Product Information:**

Description	PTFE Bottle 1000ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	1000ml
Overall Diameter:	100mm
Internal Diameter of Neck:	57mm
Overall Height:	185mm
Nominal Weight:	530g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20