

1. Part Numbers: A125-25

2. Photograph:



3. Product Information:

Description	PTFE Bottle 25ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	25ml
Overall Diameter:	33mm
Internal Diameter of Neck:	18mm
Overall Height:	61mm
Nominal Weight:	42g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 11-Apr-24