

1. Part Numbers: A125-500

2. Photograph:



3. Product Information:

Description	PTFE Bottle 500ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	500ml
Overall Diameter:	80mm
Internal Diameter of Neck:	52mm
Overall Height:	150mm
Nominal Weight:	285g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20