

Specification Number: A125-50

Revision: Date Issued:

03 27-Oct-20

- 1. Part Numbers:
 - A125-50

2. Photograph:



3. Product Information:

Description	PTFE Bottle 50ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	50ml
Overall Diameter:	43mm
Internal Diameter of Neck:	25mm
Overall Height:	76mm
Nominal Weight:	66g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Position: Date:

Kuldeep Flora **Quality Manager** 27-Oct-20



4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

SteriWare is a registered Trademark belonging to Sampling Systems Ltd.