

1. Part Number: A166-500

2. Photograph:



3. Product Information:

Description:	Powder Pot - 500ml
Material of Construction:	316L stainless steel
Nominal Volume:	500ml
Brimful Volume:	500ml
Overall Height (body & lid):	81mm
Body Height:	72mm
Maximum Body Diameter:	106mm
Lid Diameter:	120mm
Nominal Weight:	1.3kg
Finish	Bright Polish (Ra < 0.4 microns)
BSE/TSE statement:	The polishing compounds used in the manufacture of these products are of vegetable origin and so are certified BSE/TSE free

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Apr-20