



1. Part Number: A200-200

2. Photograph:



3. Product Information:

Description:	Pallet Knife with Plastic Handle
Material of Construction:	- Blade: Stainless Steel - Handle: Black PVC
Method of Construction:	Two piece construction
Surface Finish:	Better than 0.5 microns Ra
Blade Length:	200mm
Blade Width:	30mm
Overall Length:	310mm
Nominal Weight:	59g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora
Position: Quality Manager
Date: 18-Dec-24