

Specification Sheet

Specification Number: A244-310
Revision: 05
Date Issued: 30-Nov-21

1. Part Number: A244-310

2. Photograph:



3. Product Information:

Description: Industrial Liquid Funnel

Material of Construction: 304 Stainless Steel

Top Diameter: 310mm

Outer Diameter of Spout: 20mm

Overall Length: 325mm

Nominal Spout Length: 88mm

Nominal Weight: 530g

Side Handle: Spot welded to cone section of funnel

Recommended Storage Conditions: Dry and ambient temperature

Note: Due the manufacturing method there is an internal crevice where the cone and spout are joined. These funnels are primarily designed for industrial applications.

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager
Date: 30-Nov-21



4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com