

**1. Part Number:** A246-10

**2. Photo:**



### 3. Product Information:

Description	Shovel
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	1155mm
Blade Length:	350mm
Blade Width:	235mm
Nominal Weight:	2.6kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product is designed for medium duty use

### 4. Approval:

Signed:



**Name:** Kuldeep Flora  
**Position:** Quality Manager  
**Date:** 15-Feb-23