

1. Part Number: A256-05

2. Photograph:



3. Product Information:

Description:	Sampler Stands
Method of Construction:	Single piece construction
Materials of Construction:	316L stainless steel
Surface Finish:	Satin Finish (Ra < 1.0 microns)
Quantity Supplied:	2
Width:	70mm
Height:	85mm
Max. Width of 'V' notch:	40mm
Material Thickness:	2mm
Nominal Weight:	200g each
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 12-Jul-21