

## **Specification Sheet**

Specification Number: A296-4
Revision: 02
Date Issued: 2-Dec-19

**1. Part Number**: A296-4

## 2. Photograph:



## 3. Product Information:

Description Stainless Steel Bailer 4L Materials of Construction: 316L stainless steel

Nominal Volume: 4 Litres

Method of Construction: Welds ground and polished Surface Finish: Better than 0.5 microns Ra

Handle Length: 145mm

Max. Bowl Diameter: 270mm

Bowl Height: 115mm

Nominal Weight: 1010g

Graduation Marks: None

BSE/TSE statement: All polishing compounds used in the manufacture of this product are

of vegetable origin.

Recommended Storage Conditions: Dry and ambient temperature

4. Approval: Signed:

Name: Kuldeep Flora
Position: Quality Manager

Date: 2-Dec-19

