


	<h1>Product Specification</h1>	Specification Number	A324-50
		Revision	01

Part No. A324-50
Description Ladle
Nominal Volume: 50ml
Material of Construction: 304 stainless steel



Surface Finish: Better than 0.5 microns Ra
Overall Length: 275mm
Bowl Diameter: 63mm
Nominal Weight: 105g

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 22-Apr-13
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	