

1. Part Number: A326-150

2. Photograph:



3. Product Information:

Description	Ladle 150ml
Material of Construction:	316L Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	150ml
Overall Length:	370mm
Bowl Diameter:	88mm
Nominal Weight:	150g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20